

new!

Pink
pole™

infusion pump stand



infusion pump stand

Height adjustable from 70" to 120"

Weight Capacity 75 lbs. 6" off center†



P-1095-8


features

- Bag **drop prevention device** allows for controlled descent of pole.
- Offset pole sleeve welded securely to the base for reinforced stability.
- Low center of gravity design and 30 lbs. base offers improved stability.
- Heavy duty 3" sealed casters with precision ball bearings.
- Extremely durable powder coat finish base.
- Eight stainless steel IV hooks.

† 75 lbs max applied at 6" off center.
‡ 95 lbs max total pumps and bags.



Option - Handle & Tray Insert



Heavy-Duty 30 lbs. Base
26" Diameter



Offset Post

