



Cool-Touch Handle



Modular Bumper



Biohazard Indicator Wheel



Waterfall Wire Shelf

## Design

- Stainless steel construction.
- Integrated body and floor increases overall cart rigidity.
- High temp, impact resistant bumper mounted directly to cart body reducing the overall footprint.
- Cool-Touch polymer handle.
- Easy clean interior - great for hand washing and cart washing.
- 6" diameter stainless steel casters.
- CDS-342 and CDS-345 feature two direction lock casters.
- Roll-out shelves with two position latch mechanism safely locks the shelf in or out.
- Cart is designed to quickly drain during cart washing to minimize residual water.



## roboticart



**CDS-333**  
 Overall outside dimensions:  
 29½"W x 28"D x 56"H  
 Cabinet Dimensions  
 23"W x 24¾"D x 49"H  
 Usable Interior Space  
 20"W x 24¼"D x 44¼"H



**CDS-335**  
 Overall outside dimensions:  
 35"W x 32¾"D x 56"H  
 Cabinet Dimensions  
 28½"W x 29½"D x 49"H  
 Usable Interior Space  
 25½"W x 29"D x 44.¼"H



**CDS-342**  
 Overall outside dimensions:  
 46¾"W x 28"D x 39"H  
 Cabinet Dimensions  
 40¼"W x 24¾"D x 32"H  
 Usable Interior Space  
 37½"W x 24¼"D x 27"H



**CDS-345**  
 Overall outside dimensions:  
 46¾"W x 28"D x 56"H  
 Cabinet Dimensions  
 40¼"W x 24¾"D x 49"H  
 Usable Interior Space  
 37½"W x 24¼"D x 44¼"H



**CDS-356-MS**  
 Overall outside dimensions:  
 29½"W x 28"D x 39"H  
 Cabinet Dimensions  
 23"W x 25"D x 32"H  
 Usable Interior Space  
 20"W x 24¼"D x 27¼"H